

Related product : Recycling Preparative HPLC Series


 Recycling Preparative HPLC  
LaboACE LC-5060

## Recycling by GPC Column Separation of Polyvinylidene Difluoride Oligomer

### Keyword:

 Separation of PVDF, GPC Column, Size Exclusion  
Chromatography

### Introduction

In preparative HPLC, the column length is one of the key factors to get better separation. However, there is a limit in length due to restriction on the pressure the column can endure.

Recycling preparative HPLC is the solution to the problem. By cycling the sample solution back to the same column repeatedly, it causes the same effect as a longer column is used. Further, no solvent is consumed during the cycles. So it is the ideal way to efficiently increase separation (resolution) ability.

Moreover, combined use of SEC column, which separates compounds by their size, has gained great popularity among synthetic organic chemists since it can considerably save labor and time for method development as far as the sample is dissolved in some solvent.

Here is an example of recycling preparative HPLC using organic GPC column.

### Experiment & Results

Sample: Mixture of Polyvinylidene Difluoride Oligomer (n=9,10,11)

We tried to isolate the oligomers by Recycling Preparative HPLC using GPC column.

Instrument : LC-9204 (Detector : UV (254 nm))

Column : JAIGEL-2H + JAIGEL-3H in series

Mobile phase : THF

Flow rate : 3.5 mL/min




Fig. 1




Fig. 2

### Conclusion

The three oligomers were separated at the 19th cycle.