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
図7 UV trace of the triacylglycerol fraction of hydrolyzed fish oil with Lipase OF.

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POSITIVE AND NEGATIVE FAB MASS SPECTROMETRY BY FRIT-FAB LC/MS FOR THE CHARACTERIZATION OF ACYLGLYCEROLS IN FISH OIL

Masanori Hori¹, Yuko Sahashi², and Seiji Koike³

¹Nitto Technical Information Center Co., Ltd., 1-1-2, Shimozumi, Ibalaki 567, Japan

²Medical and Membrane Research Laboratory, Nitto Denko Corp., 1-1-2, Shimohozumi, Ibalaki 567, Japan

³Food, Fat and oil Laboratory, Asahi Denka Kogyo Corp., 8-4-1, Higashi-ogu, Arakawa, Tokyo 116, Japan

1. Introduction

Natural triacylglycerols (TGs), especially fish oil, are such complex mixture of different molecular species that it is difficult to separate all of them by a single chromatographic method. Therefore the experimental elucidation of the molecular species in the fish oil remains the most formidable problem in lipid chemistry. The composition analyses of TGs using FRIT-CI LC/MS with ammonia as a reagent gas have been reported previously (1), however this method was not suitable for the polar, labile and involatile fish oil. In this work FRIT-FAB LC/MS has been utilized to obtain the composition of the molecular species concerned with polyunsaturated fatty acids in the fish oil hydrolyzed with Lipase OF (*Candida cylindracea*).

2. Experimental

I

The HPLC analyses were performed on a Hewlett Packard model 1090L with reversed phase columns, ODS-Hypersil using a mobile phase of Acetone: acetonitrile (35:75) at a flow rate of 0.15 ml/min. The eluent via UV detector. Was introduced into a JEOL JMS-AX505H double focusing mass spectrometer equipped with a FRIT-FAB interface. The *m*-nitrobenzyl alcohol was used for the matrix of FAB ionization and introduced to the ion source by the post column addition.

3. Results and Discussion

According to the fundamental studies of the authentic TGs the negative ion FAB mass spectra of the TGs gave the RCOO⁻ ions as the base peak, and the other characteristic ions were not observed. Therefore the distribution of fatty acids in the oil could be determined by recording the mass chromatograms of the RCOO⁻ ions. Figure 1 shows the UV trace and the mass chromatograms of the RCOO⁻ ions of the TG fraction from the hydrolyzed fish oil (hydrolysis I-ate Was 5LL7E). The presence of eicosapentaenoic acid (EPA) and docosahexaenoic acid (DHA) in the oil was confirmed by monitoring the RCOO⁻ ions at *m/z*: 301 and 327 respectively. For example it was evident that the TG in peak 1 consisted of at least two fatty acids of EPA and DHA. The positive ion FAB mass spectrum of peak 1 (shown in Fig. 1) represents the two species of EPA-EPA-EPA (tri-EPA) and 18:4-EPA-DHA as indicated by the [MH-RCOOH]⁺ ions and the [M+H]⁺ ion. The retention time of the peak 1 also

coincided with that of the authentic tri-EPA, and the retention time of the 18:4-EPA-DHA coincided well with the value that has been calculated by the relative retention potential index theory(2) as shown in Table 1. FAB mass spectrometric techniques combined with HPLC provided detailed information on the molecular species of the TGs in the fish oil.

References

- (1) M. Hori, Y. Sahashi, and S. Koike, *Mass Spectrosc.*, 39, 133(1991).
- (2) K. Takahashi, and T. Hirano, *Nippon Suisan Gakkaishi*, 54, 523(1988).


Fig. 1. Negative ion mass chromatograms of the hydrolyzed fish oil and positive ion FAB mass spectrum of peak 1.

Table 1. TG Molecular Species of the hydrolyzed fish oil determined by FRIT-FAB LC/MS

Molecular species ^{a)}	Empirically determined retention time (min)	Predicted ^{b)} retention time (min)	Relative retention time deviation (%)
(16:4)(18:4)(EPA)	17.95	17.87	0.45
(16:4)(EPA)(DHA)	19.27	19.28	0.05
(18:4)(EPA)(EPA)	21.55	21.64	0.42
(EPA)(EPA)(EPA)	22.36	—	—
(18:4)(EPA)(DHA)	22.36	22.45	0.40
(EPA)(EPA)(DHA)	23.21	23.35	0.60
(EPA)(DHA)(DHA)	24.10	24.22	0.50

^{a)} The binding position of fatty acids residues are not discriminated.

^{b)} Values that have been calculated by relative retention potential index theory.