

## Curie Point Pyrolyzer JHP-5

### Thermal decomposition of HDPE with different Py-temp

Using setup model of curie point pyrolyzer JHP-5 with different Py temperatures.

Pyrograms are from the top of F500 (500c), F590(590c), F670(670c), and F764(764c).

From the pyrogram of F500, each peak heights are only 15% of F590, F670 and F764, polyethylene were not thermally decomposed completely.

Also, as the Py temperature goes higher to F764, it detected more small molecules.

From this analysis, optimized Py temperature is at F590 (590c).

#### Analysis Condition

Sample	HDPE (High Density Polyethylene)	Amount: about 0.3 mg	
JHP-5	Oven temperature	300°C	
	Needle temperature	300°C	
	Pyrofoil	F590	
	Pyrolysis time	5 sec	
GC-17A	Constant flow mode		
	Column flow rate	1.0 ml/min (He)	
	Velocity	36.1 cm/sec	
	Oven temperature	Initial: 40°C(3 min)	
		Rate: 10°C/min	
		Final: 320°C(69 min)	
	Head pressure	Initial: 48.9 kpa (3 min)	
		Rate: 3.8 kpa/min	
		Final: 155.3 kpa (69 min)	
	Analysis time	100 min	
	Injector temperature	320°C	
	Detector temperature	320°C	
	Equivalent time	3 min	
	Split ratio	1/100	
Total flow rate	103 ml/min		
Purge gas flow rate	5.0 ml/min (at 100 kpa)		
QP-5000	Mass range	33-550	
	EM gain	1.70 kv	
	Solvent cut	0 min	
	Scan time	0.2-100 min	
Column	DB-5MS	Inside diameter: 0.25 mm	
		Length: 30 m	
		Film thickness: 0.25 um	

# Polyethylene Thermal Decomposition Using different Py temperature

