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
Fig. 1. Chart of Peak Bands Observed in Pyrograms of 21 L- α -Amino Acids

Gly:Glycine, Ala:Alanine, Val:Valine, Leu:Leucine, Ile:Isoleucine, Phe:Phenylalanine, Tyr:Tyrosine, Trp:Tryptophan, Thr:Threonine, Cys:Cysteine, (Cys)₂:Cystine, Met:Methionine, His:Histidine, Pro:Proline, Hyp:Hydroxyproline, Asp:Aspartic acid, Asn:asparagine, Glu:Glutamic acid, Gln:Glutamine, Lys:Lysine, Arg:Arginine

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TABLE I. Degradation Products Identified in Pyrograms of Aromatic Amino Acids

No.	L-Phenylalanine	L-Tyrosine	L-Tryptophan
1	Benzene	Phenol	Indole
2	Toluene	<i>p</i> -Cresol	Skatole
3	Styrene	<i>p</i> -Hydroxystyrene	3-Ethylindole
4	Unknown	Unknown	
5	Phenethylamine		
6	Diphenyl		
7	Dibenzyl		
8	Unknown		

TABLE II. Degradation Products Identified in Pyrograms of Aliphatic Amino Acids

No.	L-Alanine	L-Leucine	L-Isoleucine
1	CH ₃ -CHO	CH ₃ -CH(CH ₃)-CH ₂ -CHO	CH ₃ -CH ₂ -CH(CH ₃)-CHO
2	CH ₃ -CN	CH ₃ -CH(CH ₃)-CH ₂ -CN	CH ₃ -CH ₂ -CH(CH ₃)-CN
3		CH ₃ -CH(CH ₃)-CH ₂ -CH ₂ -N	CH ₃ -CH ₂ -CH(CH ₃)-CH ₂ -N
		CH ₃ -CH(CH ₃)-CH ₂ -CH ₂ -N	CH ₃ -CH ₂ -CH(CH ₃)-CH ₂ -N

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
Fig.2. Pyrograms of Pepsin and Hemoglobin

TABLE III. Degradation Products Identified in Pyrograms of Pepsin

Peak No.	Degradation products	Amino acid residues assigned
1	Acetaldehyde	Alanine
2	3-Methylbutanal	Leucine
3	2-Methylbutanal	Isoleucine
4	2-Methylpropanenitrile	Valine
5	Toluene	Phenylalanine
6	2-Methylbutanenitrile	Isoleucine
7	3-Methylbutanenitrile	Leucine
8	Pyrrole	Proline, hydroxyproline
9	Phenol	Tyrosine
10	<i>p</i> -Cresol	Tyrosine
11	Unknown	Glutamic acid, glutamine
12	Unknown	Aspartic acid, asparagine
13	Indole	Tryptophan
14	Skatole	Tryptophan

a) Peak numbers correspond to those in the pyrogram of pepsin


Fig.3. Pyrograms of Wool Fiber, Silk Fiber and Leather

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
Fig.4. Pyrogram of Wool Fiber (Raw fiber,KBr disk of raw fiber, Ash ball of fiber and KBr disk of ash ball)


Fig.5. Plots of N(Pro),N(Tyr) and N(Trp) on PGC of Wool, Silk,Leather and Nylon 6,6

○:raw fiber, ●:KBr disk of raw fiber,
 ▲:ash ball of fiber, ▲:KBr disk of ash ball,
 □:N value determined using reported values

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
Fig.6. Analytical Procedure for Quaternary Ammonium Compounds in Human Urine


Fig.7. Total Ion Collector of Eight Quaternary Ammonium Compounds

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
Fig.8. Pyrograms of methylbenactyzium bromide

- (a) urine sample effluent, directly pyrolyzed
- (b) urine sample effluent, pyrolyzed after PTLC
- (c) urine control effluent, pyrolyzed after PTLC

Table IV. Comparison of PGC and PTLC-PGC for 7 Quaternary Ammonium Compounds

Compounds	PTLC-PGC/PGC (%)
Distigmin bromide	150
Bethanechol chloride	n.d.
Pancuronium bromide	807
Propantheline bromide	341
Suxamethonium chloride	65
Neostigmine methylsulphate	76
Benzethonium chloride	25

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